

MX2

PLUS



THE LITTLE PRODIGY.

The use of the latest generation magnets enhance the performance of the MX2 Plus micromotor.

The result is unparalleled balance and comfort, more power and perfect control on speed and torque.

**« GREATER POWER,
MORE TORQUE. »**



POWER UNDER CONTROL.

A major attribute of Bien-Air micromotors is their elevated power output with constant torque for effective and effortless use.



infinitely
small.

**12%
LIGHTER***



**15%
SHORTER***

*Than MX2



« A NEW OPTIMIZED DESIGN. »

**10.000
WORKING CYCLES.**

MX2 PLUS has been successfully tested at 10.000 working cycles, making proof of withstanding repeated sterilization cycles to respect the strictest hygiene standards.



ACCURATE CONTROL.

The proprietary Smart Logic™ electronic control provides Bien-Air micromotors a step-less and vibration-less rotation, allowing for the most meticulous and precise dental work.

LED LIGHT.

Optimized light intensity and colour temperature for better comfort during your operations.

MAINTENANCE-FREE.

Bien-Air micromotors' ball bearings are greased for life, making no further greasing necessary.



MX2 PLUS

- Dimensions
Ø 21 x 44 mm
- LED with adjustable intensity
- Lifetime lubricated ball bearings
- Stainless steel
- 3-year warranty

MX2 PLUS REF 1600677-001

WWW.BIENAIR.COM

S W I S S  M A D E

More information

